

γstream

Active, stand-alone, fully featured
MCA tube base
for scintillation
spectroscopy



Gamma Spectroscopy goes wireless with the innovative CAEN tube base digital MCA: Untie your detector and get on-line!

Features

- Fully stand-alone MCA including high voltage power supply, preamplifier, battery and data-storage on SSD.
- Compatible with scintillation detectors as NaI(Tl) and LaBr₃(Ce) using standard 14-pin and 10-stage PMTs
- Different acquisition modes available: PHA, PHA with time stamp, MCS, SCA and signal inspector
- Software selectable coarse and fine gains
- ARM based CPU for unattended operations
- Gain stabilizer based on natural or calibration radioactivity
- Wire and wireless connectivity through USB, Ethernet, Bluetooth and WiFi
- Front panel cover against dust and rain for outdoor operation
- Front panel auxiliary digital I/O connectors for synchronization, external trigger, coincidence/ anticoincidence modes, veto, MCS, SCA
- Supported by CAEN MC²Analyzer software GUI
- Supported by GammaTOUCH application for smartphones or tablets with Android OSs

γstream is an active MCA integrated in a 14-pin PMT base for gamma-ray spectroscopy with scintillation detectors, suited for portable and mobile spectrometry. It integrates high voltage power supply, preamplifier and digital MCA. The embedded CPU runs an OS able to execute custom routines for automatic wireless spectroscopic analysis or data monitoring.

Featured with battery and data-storage it is able to run without external connection or control. The flexibility on data logging modes makes *γstream* suitable for most survey situations.

It may acquire and record data in different modes: PHA, time stamped list mode, MCS, SCA, signal inspector. *γstream* can also operate outdoor thanks to a front panel cover protecting from water and dust.

Wireless and wired connection interfaces are available: WiFi, Bluetooth, USB and Ethernet. *γstream* can also be integrated in more complex systems thanks to auxiliary I/O connectors.

γstream is suited for a variety of environment thanks also to the software suites GammaTOUCH and MC²Analyzer, both provided with user-friendly GUIs. GammaTOUCH can run on smartphones or tablets with Android OSs, providing an immediate and easy data acquisition control. It provides georeferenced and time-stamped histograms logging the GPS location and time of the mobile device. On the other side MC²Analyzer can be run on recent Windows, providing the user with histogram analysis tools.

Software

GammaTOUCH


GammaTOUCH is a software application for Android OSs running on personal devices like smartphones or tablets. It is designed as an easy-to-use assistant for outdoor operations and provides an immediate control of *γstream*. Energy spectra can be displayed, and statistics can be monitored runtime. Through the GammaTOUCH it is possible to integrate the location and time of the mobile device GPS and the time-stamped histogram.



γstream rear view



All CAEN Control Software are available for **free download** on the web site.



Typical NaI detector with 14-pin PMT. *γstream* is plugged on the PMT tube.



News from Catalog web page
www.caen.it/news



n Small details
Great differences

n CAEN
Tools for Discovery

Copyright © CAEN SpA - 2015
All rights reserved. Information in this publication supersedes all earlier versions. Specifications subject to change without notice.
Printed in March 2015 - ADOCUME00112 - BF3264 - rev03

CAEN SpA
Via Vetraia 11
55049 - Viareggio - Italy
Phone +39.0584.388.398
Fax +39.0584.388.959
info@caen.it
www.caen.it

CAEN GmbH
Klingenstraße 108
42651 - Solingen - Germany
Phone +49.212.2544077
Fax +49.212.2544079
info@caen-de.com
www.caen-de.com

CAEN Technologies, Inc.
1140 Bay Street - Suite 2C
Staten Island, NY 10305 - USA
Phone +1.718.981.0401
Fax +1.718.556.9185
info@caentechnologies.com
www.caentechnologies.com